Z_M1

Conical implants with external hex connection





Z_M1

Conical implants with external hex connection





Important information

Please read carefully before using Ziacom® products

General information

This document contains basic information on the use of original Ziacom® dental implant systems, hereafter referred to as "Ziacom® dental implants" or simply "Ziacom® products". This document has been created as a quick guide for clinicians responsible for treatment, hereafter the "user", and therefore, is neither an alternative nor a substitute for specialised training or professional clinical experience.

Ziacom® products must be used according to a suitable treatment plan and in strict compliance with the manufacturer's surgical and prosthetic protocols. Carefully read the product-specific surgical and prosthetic protocols and the instructions for use and maintenance before using any Ziacom® product. You can find this information on our website, www.ziacom.com, or request it from your nearest authorised Ziacom® distributor.

Liability, safety and warranty

The instructions for the use and handling of Ziacom® products are based on internationally published literature, current clinical standards and our clinical experience so they should be understood as general guidance. The handling and use of Ziacom® products is the sole responsibility of the user as it is outside the control of Ziacom Medical SL. Ziacom Medical SL, its subsidiaries and/or its authorised distributors disclaim all responsibility, whether explicit or implicit, total or partial, for possible damage or injury caused by poor handling of the product or any other situation not considered in their protocols and manuals for the correct use of their products.

The user must ensure that the Ziacom® product is appropriate for the intended procedure and end purpose. Neither these instructions for use nor the work or handling protocols for the products release the user from this obligation. Ziacom® products must be used, handled and applied by clinicians with the appropriate training and qualifications required according to current legislation in each country.

The total or partial use, handling and/or application of Ziacom® products at any stage of their implementation by personnel who are unqualified or lack the necessary training will automatically void any type of warranty and may cause severe damage to the patient's health.

Ziacom® products are part of their own system, with their own design characteristics and work protocols, including dental implants, abutments or prosthetic components and surgical or prosthetic instruments. The use of Ziacom® products in combination with elements or components from other manufacturers could result in treatment failure, damage to tissues or bone structures, inadequate aesthetic outcomes and severe damage to the patient's health. Therefore, only original Ziacom® products should be used.

The clinician in charge of the treatment is solely responsible for ensuring the use of original Ziacom® products and that they are used according to the corresponding instructions for use and handling protocols throughout the implant procedure. The use of any other non-original Ziacom® components, instruments or products, whether alone or in combination with any original Ziacom® products, will immediately void the warranty of the original Ziacom® products.

See the Ziacom Medical SL. Warranty Programme (available on the website or by contacting Ziacom Medical SL, its subsidiaries or authorised distributors).

Warning. Not all Ziacom® products are available in all counties. Check availability in your country.

The Ziacom® brand and the names of other products and services, including their logos, that are mentioned in this document or on the website www.ziacom.com are registered trademarks of Ziacom Medical SL.

Ziacom Medical SL reserves the right to modify, change, remove or update any of the products, prices or technical specifications referenced on this website or in any of its documents without prior notification. All rights reserved. The reproduction of this document, whole or in part and in any medium or format, without the corresponding written authorisation from Ziacom Medical SL is prohibited.



Together for health




Table of contents

Together for health

Ziacom® quality

Zitium® Titanium

06

06

06

06

Investment in innovation and training	07
Ziacom® across the globe	07
- Regional headquarters	07
- Subsidiaries	07
ZM1 conical implants with external hex connection	
ZM1 implant	10
Characteristics	10
Diameters and lengths	11
Surface treatments	12
- Titansure surface treatment	12
Product presentation	14
ZM1 references	16
Recommendations for use	17
How to use this catalogue	18
Product sheet	18
Symbology	18
Abutments Direct-to-implant restorations	20
Abutments Transepithelial restorations	30
Surgical instruments	40
Prosthetic instruments	48
Surgical protocol	54
Cleaning disinfection and sterilisation	74

The Company

Together for health

Ziacom® has been working for more than 20 years to improve the oral health and well-being of patients around the world by designing and manufacturing innovative, high-quality dental implant, prosthetic component, surgical instrument and biomaterial solutions.

The company was founded in 2004 with 100% Spanish capital and began its activity as a manufacturer of dental implants and attachments for several European companies before later launching its own **brand of implant systems** in 2006.

In 2015. Ziacom® introduced its diversification strategy with the development of **new business lines** and new product lines and the launch of a **new portfolio**, which helped the company achieve a **15%** share of the Spanish market in 2016 with the sale of more than 230.000 implants.

In 2022, the company began an ambitious growth plan with new goals of international expansion, broadening and diversification of its portfolio of products and services and a change in corporate identity.

Ziacom® quality

Commitment to quality and innovation has been part of the values and the essence of Ziacom® since the beginning.

That is why we use state-of-the-art technology in every stage of our products' production cycle, from design and manufacture to quality assurance, cleaning and packaging. All of our products are also manufactured using only high-quality raw materials after applying strict controls to select our main suppliers.

Ziacom Medical S.L. is a licensed manufacturer of medical devices and an AEMPS (Spanish Agency of Medicines and Medical Devices) 6425-PS marketing authorisation holder. Our quality management system is certified in accordance with the requirements of ISO standards 9001:2015 and 13485:2018. and is also GMP 21 CFR 820 compliant.





Thanks to our ceaseless endeavours to offer our clients unsurpassable quality, all our implants have a lifetime guarantee.


See the General Conditions for Accessing the Warranty for Ziacom® products.

7itium® titanium

Ziacom® ZM1 implants are made from extra-high-strength grade 4 Zitium® titanium, which bestows them with substantially improved elastic limit and mechanical properties.

With **Zitium**® our implants meet the requirements of standards ASTM F67 and ISO 5832-3. and are certified in accordance with EU Regulation 2017/745. attaining the corresponding CE marking from notified body 0051.

Properties of Zitium® titanium















Ziacom® dental implants are all sterilised using beta-ray radiation at 25 kGy, apart from the DSQ orthodontic implants, which are supplied non-sterile.

All the products (except dental implants) listed in this Ziacom® catalogue are supplied non-sterile and must be sterilised



i



Investment in innovation and training

In order to always offer the very best solutions for the **well-being of every patient**, and thanks to the experience and dedication of our **highly-qualified professionals** and **innovative Technological Centre**, our R&D&I team works incessantly in the field of **research and innovation** to **improve** our products and develop **new solutions** to meet the demands and needs of both patients and dentists.

We also invest in **research** and **ongoing training** as a way of providing **scientific support to the sector** and we firmly believe in training **young professionals** to best ensure **advances in the dentistry field**.

We therefore work closely with **training centres**, **universities and** scientific bodies to create a practical and specialised teaching

environment to promote and strengthen their knowledge, abilities and professional growth.

In order to enhance our investment in the training and **development of dental professionals**, we have **specific areas at our facilities for hands-on training** and **practicals**, **state-of-the-art** training equipment and also a **physical and virtual showroom** where professionals can see all our dental solutions first hand.

Ziacom® across the globe

We are committed to making oral health available to patients all over the world and have a solid **internal growth and expansion plan** to increase the company's **international presence** in those **areas where our products are already well-established** and to **expand into new areas**

In order to achieve this, we offer our **international associates** a **trusting and collaborative** partnership by adapting to their **local needs** and providing solutions that are specific to each market.

As part of our commitment to meet the specific **quality, regulatory and legal requirements of each country**, for both the registration and distribution of our products, we have **specific certifications** from each of the countries in which we trade.

Regional headquarters

Ziacom Medical

Madrid - SPAIN Calle Búhos, 2 - 28320 Pinto ♣34 91723 33 06

Subsidiaries

Ziacom Lusobionic

Av. Miguel Bombarda, 36 - 5° B 1050 -165 - Lisbon - PORTUGAL \$\tilde{\text{\$\chi}}\ +351 215 850 209

Ziacom Medical USA LLC

Miami - USA 333 S.E 2nd Avenue, Suite 2000 Miami, FL 33131 - USA

+1(786) 224 - 0089


Ziacom ITS

Viale del Lavoro, 14 35010 Vigonza Padova - ITALY \$\infty\$ +39 049 603310

Ziacom Safe implant

Av. Iñaquito, Edificio Metropolitan, Oficina 304 170507, Quito - ECUADOR \$\infty\$ +593 96 368 0879

Please see the up-to-date list of Ziacom® distributors at www.ziacom.com or email us at export@ziacom.com



ZM1

ZM1 Conical implants with external hex connection



ZM1 implant

Characteristics

CONNECTION

- External hex connection: simple and versatile.
- Screw channel with upper guide: facilitates screw insertion.

CORTICAL ZONE


- · Microthread design: preserves marginal bone.
- · Microthread extension: improves load distribution.
- · Macrodesign: optimal cortical compression.

BODY


- · Reduced-angle active threads: improve stability during insertion and increase BIC (bone-to-implant contact).
- · Double threaded: quick insertion and shorter surgical time.
- Self-tapping active apex: facilitates insertion with underdrilling.
- Transverse apical windows: collect remnants of bone during insertion.
- · Optimised morphology: high primary stability.
- · Atraumatic apex: no damage to anatomical structures.

CONICAL DESIGN

- · Facilitates shaping in low density bone.
- · Indicated for immediate loading.
- Indicated for cases of apical convergence and/or collapse.



Dimensions of the implant collar



Advantages

- · The external connection has been used for many years by implantologists, with an acceptable degree of confidence.
- · Having been on the market for many years, there is a large number of studies to back up its efficacy.
- · It is a good option when there is some divergence between implants.
- · Suitable for multi-unit restorations.
- · Very easy to restore.
- · It has a wide range of options for restoration.

Ziacom®



Diameters and lengths

					LENGTH (L)			
Ø DIAMETER	Ø PLATFORM	6	7	8.5	10	11.5	13	14.5
NP 3.30	3.30							
RP 3.60								
RP 4.00	4.10							
RP 4.40								
WP 4.80	5.00							

ZM1

Dimensions in mm.

ZM1 implant

Surface treatments

■ Titansure surface

Implants inserted following surface treatment are known to benefit from improved osseointegration by increasing the bone-to-implant contact area. This is partly due to the implant's chemical composition and topographical characteristics.


With its **Titansure** surface treatment, Ziacom® achieves contaminant-free surface topography and optimal average macro and microporosity values, which are key specifications for achieving prompt and proper osseointegration and, in turn, extremely reliable and predictable implants.


■ ANALYSIS OF THE TITANSURE SURFACE TREATMENT


Titansure is an SLA surface treatment created through a subtraction process involving sandblasting with white aluminium oxide and double acid-etching with hydrofluoric acid and a sulphuric/phosphoric acid mix.

Surface morphology analysis

With the aid of a scanning electron microscope (FEI TENEO, Thermo Fisher Scientific Inc., Waltham, MA, USA), we can see the rough, porous surface creating numerous cavities with thin, sharp edges.








Surface elemental analysis

We used an energy-dispersive X-ray spectrometer (Octane Super, Edax-Ametek, Mahwah, NJ, USA) to analyse the chemical composition at the surface.



Compositional analysis of implant surface

ELEMENT	WEIGHT (%)
CK	9.32 (10.23)
AL K	-
Ti K	89.53 (11.77)

No aluminium was detected

Results are expressed as the mean and standard deviation of the mass percentage (WEIGHT (%)).




Surface roughness analysis

The roughness study was conducted with a Sensofar S NEOX interferometric-confocal microscope (Sensofar Medical, Terrasa, Spain) and SensoMAP Premium 7.4 software. The quantitative roughness profile parameters applied were: average roughness (Ra), root-mean-square roughness (Rq), maximum profile peak height roughness (Rp) and maximum profile valley depth roughness (Rv).

Ra (µm) (SD)	Rq (µm) (SD)	Rp (µm) (SD)	Rv (µm) (SD)
0.82 (0.10)	0.97 (0.08)	1.84 (0.04)	2.21 (0.01)

The 3D surface roughness (Sa), 3D root mean square height (Sq), maximum 3D peak height (Sp) and maximum 3D pit depth of the selected area (Sv) were also recorded.

Sa (µm) (SD)	Sq (µm) (SD)	Sp (µm) (SD)	Sv (µm) (SD)
0.76 (0.01)	0.97 (0.01)	4.20 (0.12)	4.62 (0.20)



The data were extracted from:

Rizo-Gorrita, M.; Fernandez-Asian, I.; Garcia-de-Frenza, A.; Vazquez-Pachon, C.; Serrera-Figallo, M.; Torres-Lagares, D.; Gutierrez-Perez, J. Influence of Three Dental Implant Surfaces on Cell Viability and Bone Behavior. An In Vitro and a Histometric Study in a Rabbit Model. Appl. Sci. 2020, 10(14), 4790

OPTIMAL OSSEOINTEGRATION

The **Tibansure** surface has a three-dimensional surface structure with high peaks and broad troughs, which is known to be highly effective at promoting the coagulation cascade and the release of growth factors through platelet activation [Kim, H.; Choi, S.H.; Ryu, J.J.; Koh, S.Y.; Park, J.H.; Lee, I.S. The biocompatibility of SLA-treated titanium implants. Biomed. Mater. 2008, 3, 025011.].

This type of surface may have an osteogenic effect thanks to its different topographical features at a micrometer and nanometer level, which has a very similar morphology to the osteoclastic bone resorption cavities [Le Guehennec, L.; Goyenvalle, E.; Lopez-Heredia, M.A.; Weiss, P.; Amouriq, Y.; Layrolle, P. Histomorphometric analysis of the osseointegration of four different implant surfaces in the femoral epiphyses of rabbits. Clin. Oral Implants Res. 2008, 19, 1103–1110].

For more information on the surface treatment, please see the literature available atwww.ziacom.es/biblioteca



ZM1 13 Z

ZM1 implant

Product presentation

Blister packaging

Available for implants with **Tibansure** surface. The blisters are heat-sealed and include identification labels for product traceability. There is a flap for easy opening in the surgery while preventing accidental opening.




IMPORTANT

Do not open the sterile container until just before inserting the implant.

Outer identification label

Ziacom® implants are supplied in a sealed cardboard box that includes a product identification label with a description of their main characteristics.



Description of the symbology used

C ∈ § CE marking (MDR) and notified body number. Do not use if package is damaged. MD Medical device symbol. Single-use product. # Model code. Consult instructions for use. REF Product name. Product use-by date. LOT Product batch number. Date of manufacture. Unique device identifier STERBLE R Sterilised by radiation. TT Titansure surface treatment. X Temperature limit. Titansure Active surface treatment. Caution, consult attached documentation. Rx Only Prescription only. Do not resterilise. Product distributor.

For full details on the product presentation and instructions for use (IFU), go to www.ziacomes/ifus or scan the QR code on the box.



Z 14 Ziacom®



■ ZPlus Mount

The ZM1 implant includes the **ZPlus** mount, a multi-functional abutment made in grade 5 ELI titanium (medical grade), which allows easy handling of the implant during the surgical procedure and incorporates several usage functions. Additionally, the concept of the **ZPlus** Mount is based on reducing treatment costs, as it works equally well as an implant mount, impression abutment, or abutment for provisional cemented or screwed restorations.

The ZPlus Mount is available in the Zinic® SX, Zinic® MTX, ZM4, ZM4 MT and ZM1 ranges.

As indicated, the **ZPlus** Mount can be used as a provisional abutment. In such cases, the **ZPlus** should be sculpted extra-orally and adjusted on an analogue – preferably a lab model or clamp. Check also the structural integrity of the mount and screw, to ensure that they have not suffered any deformation or damage due to excessive insertion torque or forced removal manoeuvre. Additionally, verify on an analogue that the **ZPlus** fixing screw is well fitted and that the connection is secure.

IMPORTANT


Always follow the surgical protocol when placing the implant. This will protect the mount and screw from possible damage which could prevent it being used later as an impression or provisional abutment. Use each **ZPlus** only with the implant to which it belongs. To avoid mix-ups, keep the **ZPlus** and screw with the patient's ID, listing the corresponding reference and batch number. The **ZPlus** has 3 flat sides. After finishing the implant placement procedure, ensure that one of these faces into the vestibular cavity.









As an impression transfer



As a provisional abutment for cemented or screwed restorations





ZM1



15 7


ZM1 implant

ZM1 references

■ Specifications of ZM1 with ZPlus - Titansure

IMPLANT

	Ø Size (mm)	Ø Core (mm)	Length (mm)	Ref. Titansure	
			8.5	ZM13385	,610.
			10.0	ZM13310	
3.30	3.30	2.80/1.70	11.5	ZM13311	
<u>:</u>	2.00		13.0	ZM13313	晉
١			14.5	ZM13314	晋
			8.5	ZM13685	All lines
	3.60 3.10/1.80	10.0	ZM13610		
		11.5	ZM13611		
		13.0	ZM13613	冒	
			14.5	ZM13614	=
			6.0	ZM14006	
			7.0	ZM14007	
		4.00 3.40/2.10	8.5	ZM14085	
	4.00		10.0	ZM14010	
			11.5	ZM14011	
		13.0	ZM14013	臺	
			14.5	ZM14014	
			6.0	ZM14406	
			7.0	ZM14407	
			8.5	ZM14485	
	4.40	3.80/2.30	10.0	ZM14410	
			11.5	ZM14411	冒
			13.0	ZM14413	=
			14.5	ZM14414	
			6.0	ZM14806	.000
			7.0	ZM14807	
	4.80	4.10/2.40	8.5	ZM14885	
	4.00	7.10/2.70	10.0	ZM14810	
			11.5	ZM14811	臺
			13.0	ZM14813	






^{*} Screw included with each implant.


Metric




Metrics 1.80 (NP) and 2.00 (RP/WP).

Platform







(1) External hex height. (2) Distance between faces of the external hex. (3) Diameter of working platform.



Recommendations for use

All implant treatments must respect the natural biomechanical stability of the oral cavity and allow the natural emergence of the dental crown through the soft tissue. The implantologist must assess the quantity and quality of bone currently in the implant area and consider the need for prior or simultaneous bone regeneration, as appropriate.

Ziacom® has a wide range of implants available to cover every restoration possibility. The squares on the odontogram shown represent the implant diameters and platforms recommended for each tooth position.

These recommendations are valid for the replacement of teeth with single restorations, bridges, hybrid work or overdentures.

Remember to maintain minimum distances between adjacent implants and between implants and teeth in order to preserve interdental papilla, bone vascularisation and natural emergence profiles.

Selection of the appropriate implant for each case is the sole responsibility of the implantologist. Ziacom® advises all clinicians to take into account the warnings based on scientific evidence which can be found in the product catalogues and our website.

■ CLARIFICATIONS ON DRILLING MEASUREMENTS AND TECHNIQUES

- IMPLANT SIZE: identifies the diameter and length of the implant.
- IMPLANT BODY: diameter of the implant core.
- DRILL SIZE: diameter and length of the drill bit.
- DRILLING TECHNIQUE: We have developed various drilling protocols as a blueprint for dealing with different situations that arise when performing implant surgery.

For more information on implant size selection, see the literature available at www.ziacom.com/biblioteca



Dental chart

Z_M1

Implant diameter(1)

NP

A RP

B RP

c RP

Ø3.30 mm Ø3.60 mm Ø4.00 mm Ø4.40 mm. Ø4.80 mm (1) Diameters available for analogue platforms.


Implant crown diameter

NP

RP

WP

Ø3.30 mm Ø4.10 mm Ø5.00 mm




IMPORTANT

Short, 6.00 mm and 7.00 mm implants are ONLY recommended for use in combination with normal length splinted implants (≥ 10.00 mm).

ZM1

How to use this catalogue

Product data sheet



Symbology

Symbol	Meaning	Symbol	Meaning	Symbol	Meaning
ROT	Rotatory element		Tx30 connection	Steel	Made from steel
NO	Non-rotatory element	MX,XX	Size in millimetres	Co-Cr +castable	Made from cobalt-chromium + castable plastic
	Use with manual torque (see table on p. 38)	45°	45° screw support	Cobalt Chromium	Made from cobalt-chromium
XX	Maximum operating torque	90°	90° screw support	PEEK	Made from PEEK
Ncm 10 20 30 40 50 60 70	Ratchet torque range		Use in rotation with a CA	Full	Made from castable plastic
Galaxy	Galaxy connection	XX Rpm	Maximum rotation speed	Plastic	Made from plastic
1,25mm	Screw connection	XX USES	Maximum number of uses	\$\$\$	Recommended sterilisation temperature
Kirator	Kirator connection		Single-use product	Non sterile	Unsterilised product
Nature	Nature connection	Grade 5 ELI Titanium	Made from grade 5 ELI (extra-low interstitial) titanium		Use with abundant irrigation
Basic	Basic connection	Grade 2 Titanium	Made from grade 2 titanium	∑XX _○	Maximum angle
XDrive	XDrive connection	Stainless Steel	Made from stainless steel		

Abutments Direct-to-implant reconstructions



Direct-to-implant restorations

2nd STAGE AND IMPRESSIONS

Healing abutment





Platf.	Height (H)	Reference
	1.50	HAEX2015
	3.00	HAEX2030
	5.00	HAEX2050
	7.00	HAEX2070
	1.50	HAEX3415
	3.00	HAEX3430
	5.00	HAEX3450
	7.00	HAEX3470
	1.50	HAEX5015
	3.00	HAEX5030
	5.00	HAEX5050
	7.00	HAEX5070

Anodised NP RP WP













Impression abutment





Platf.	Height (H)	Reference
	11.80	TCEX2011
	11.80	TCEX3411
	7.80/Short	TCEX3402
	11.80	TCEX5011
	7.80/Short	TCEX5002

Anodised NP RP WP



Impression abutment screw





Platf.	Height (H)	Reference
	0.00	LTSS2000
	3.00	LTSS2001
	6.00	LTSS2002
	9.00	LTSS2010
	0.00	LTSS3400
	3.00	LTSS3401
	6.00	LTSS3402
	9.00	LTSS3410
	0.00	STSS3400*

Anodised NP RP/WP



*Screws to take fast impressions with short impression transfer.

Impression abutment screw - Quickly Screws



Platt.	Height (H)	Reference
	3.00	LT2001
	6.00	LT2002
	3.00	LT3401
	6.00	LT3402


Anodised NP RP/WP



Height (H) is calculated with regard to normal impression abutment height. When using a short impression abutment, consider the difference between the heights of the abut-








Anodised |



1.60

Pick-up impression transfer



Pack of 4 units. DO NOT sterilise in an autoclave. Drillable.

Z2Plus Snap-On impression abutment





PUEX5000

Platf.	Height (H)	Reference
	1.50	Z2NPEX10
	1.50	Z2RPEX10
	1.50	Z2WPEX10

Anodised NP RP WP



IMPORTANT
Use the laboratory screw to tighten this impression abutment.

Z2Plus Snap-On impression transfer





Platf.	Height (H)	Reference
	8.00	ZPU3400
	8.00	ZPU5000



Pack of 4 units. DO NOT sterilise in an autoclave. Drillable.

Implant analogue





Platf.	Length (L)	Reference
	12.00	IAEX2000
	12.30	IAEX3400
	12.00	IAEX5000



3D implant analogue - Individual





Platf.	Length (L)	Reference
	13.00	IAEX2008D
	13.00	IAEX3408D
	13.00	IAEX5008D



3D implant analogue - Pack



	7 7	
Platf.	Length (L)	Reference
	13.00	IAEX2008DC*
	13.00	IAEX3408DC*
	13.00	IAEX5008DC*
NO L1,25mm M1,80 M2,00 Stainless Steel		

^{*}Includes base screw Ref. DSIADI and lateral screw Ref. DSIADT for analogue connection.

Screws - 3D analogue





Type	Length (L)	Reference
Base screw (1)*	-	DSIADI
Lateral screw (2)*	15.00	DSIADT




*Pack of 4 units.

21 🗷 ZM1

FIXING ELEMENTS

Clinical screw





Kiran clinical screw

For ZiaCam Ti-Base or metal structures

Length (L)

8.30

8.30

M1,80 M2,00

Special Kiran screw with surface treatment



Platf.



Reference

DSEI2010

DSEI3410



Laboratory screw



Platf.	Length (L)	Reference
	7.40	LB102000
	7.40	LB103400

NOT suitable for use as the final clinical screw.

Kiran Tx30 clinical screw





For ZiaCam Tx30 abutments and Ti-Base

Platf.	Length (L)	Reference
	6.80	DSEI2010TX
	6.80	DSEI3410TX











Special Kiran screw with surface treatment Use only with Tx30 screwdrivers

PROVISIONAL

Provisional abutment





Rotatory

Platf.	Length (L)	Reference
	9.50	RUEXT2010
	9.50	RUEXT3410
	9.50	RUEXT5010

Anodised NP RP WP



Non-rotatory

Platf.	Length (L)	Reference
	9.50	NUEXT2010
	9.50	NUEXT3410
	950	NUFXT5010

Anodised NP RP WP



Provisional abutment

Aesthetic and immediate loading abutments





Rotatory

Platf.	Length (L)	Reference
	9.50	RUEXP2010
	9.50	RUEXP3410
	9.50	RUEXP5010



Non-rotatory

Platf.	Length (L)	Reference
	9.50	NUEXP2010
	9.50	NUEXP3410
	9.50	NUEXP5010

Ziacom®





SCREWED UCLA ■ MACHINED BASE UCLA Machined base abutment UCLA + Castable abutment Rotatory Rotatory Platf. Length (L) Reference 11.00 RUEX2000 Platf. Length (L) Reference 11.00 RUEX3400 10.60 BRUEX20 11.00 RUEX5000 10.60 BRUEX34 10.60 BRUEX50 ROT ROT Non-rotatory Platf. Non-rotatory Length (L) Reference 11.00 NUEX2000 Platf. Length (L) Reference 11.00 NUEX3400 10.60 BNUEX20 11.00 NUEX5000 10.60 BNUEX34 10.60 BNUEX50

ZM1 23 **Z**

SCREWED

■ Tx30 VARIABLE ROTATION ABUTMENT

Mach. base abutment Tx30

+ 2 castable abutments (15° and 20°)



Mach, base abutment Tx30

+ 2 castable abutments (20° and 25°)



Rotatory

Platf.	15° Length (L)	20° Length (L)	Reference
	11.40	11.20	BRUEX20TX
	11.40	11.20	BRUEX34TX
	11.40	11.20	BRUEX50TX





Rotatory

Platf.	20° Length (L)	25° Length (L)	Reference
	11.20	11.00	BRUEX20TX1
	11.20	11.00	BRUEX34TX1
	11.20	11.00	BRUEX50TX1





Non-rotatory

	Platf.	15° Length (L)	20° Length (L)	Reference
		11.40	11.20	BNUEX20TX
		11.40	11.20	BNUEX34TX
		11.40	11.20	BNUEX50TX





Non-rotatory

Platf.	20° Length (L)	25° Length (L)	Reference
	11.20	11.00	BNUEX20TX1
	11.20	11.00	BNUEX34TX1
	11.20	11.00	BNUEX50TX1
NO 30 M1,80 M2,00 45 Co-Cr 4cstable			



Includes special Kiran Tx30 screw with surface treatment Ref. DSEI2010TX (NP)/DSEI3410TX (RP/WP)) for all Tx30 Variable Rotation abutments.

■ TX30 VARIABLE ROTATION ABUTMENT

The Tx30 variable rotation abutment comprises a Cr-Co machined base that accepts 15°, 20° or 25° angled castable abutments and a Kiran clinical screw with a special Tx30 connection.

The Cr-Co base ensures a perfect fit and seal with the implant connection and the different angles of the castable abutments can be used to choose the best position for the correct emergence of the restoration screw access channel.



Identifying grooves for the castable angles





CEMENTED

Straight Abutment





Straight Abutment





Platf.	Height (H)	Reference
	1.50	STAEX2015
	2.50	STAEX2025
	3.50	STAEX2035
	1.50	STAEX3415
	2.50	STAEX3425
	3.50	STAEX3435
	1.50	STAEX5015
	2.50	STAEX5025
	3.50	STAEX5035

Platf.	Height (H)	Reference
	1.50	STEX2015
	2.50	STEX2025
	3.50	STEX2035
	1.50	STEX3415
	2.50	STEX3425
	3.50	STEX3435
	1.50	STEX5015
	2.50	STEX5025
	3.50	STEX5035











Anodised NP RP WP



15° angled abutment





25° angled abutment





Platf.	Height (H)	Reference
	1.50	A1EX2015
	2.50	A2EX2015
	1.50	A1EX3415
	2.50	A2EX3415
	1.50	A1EX5015
	2.50	A2EX5015

Plati.	Height (H)	Reference
	1.50	A1EX2025
	2.50	A2EX2025
	1.50	A1EX3425
	2.50	A2EX3425
	1.50	A1EX5025
	2.50	A2EX5025











25 2 ZM1

Direct-to-implant restorations

OVERDENTURES

Kirator



Kirator abutment

Platf.	Height (H)	Reference
	1.00	LOEX2001
	2.00	LOEX2002
	3.00	LOEX2003
	4.00	LOEX2004
	5.00	LOEX2005
	6.00	LOEX2006
	1.00	LOEX3401
	2.00	L0EX3402
	3.00	LOEX3403
	4.00	LOEX3404
	5.00	LOEX3405
	6.00	LOEX3406
	1.00	LOEX5001
	2.00	LOEX5002
	3.00	LOEX5003
	4.00	LOEX5004

Golden surface treatment.













Includes the Kirator abutment with sterilisable polyoxymethylene applicator (Tecaform AH-POM-C).

Related abutments

Kirator impression transfer







Kirator

Kirator analogue



Reference IATORK01

System	Height (H)	Reference
Kirator	6.50	TCRK3400





Pack of 4 units. DO NOT sterilise in an autoclave. Drillable.

Kirator processing kit



plastic caps.



Reference

-	

and pink plastic caps.

2.05
Titanium housing

Reference

Kirator processing kit	TP8520
irator processing kit consisting of: Titaniu	

System

Sterilise the metal coping in an autoclave. Plastic caps and the disc should be cold sterilised. See Instructions for Cleaning and Sterilising on the Ziacom® website.

System	Retention (kg)	Reference
	Light/1.20 kg	TPK100
Kirator	Standard/1.80 kg	TPK200
	Strong/2.70 kg	TPK300

Pack of 4 plastic Kirator retainer caps.



NOT autoclavable – use cold steriliser. Maximum divergence of 22° between implants.

Kirator divergence processing kit

Length (L)

13.00



System

3,300	r terer er ice
Kirator processing kit	TP8520D
Kirator divergence processing kit comprisi	

Sterilise the metal coping in an autoclave. Plastic caps and the disc should be cold sterilised. See Instructions for Cleaning and Sterilising on the Ziacom® website.

System	Retention (kg)	Reference
	Light/1.20 kg	TPK110
Kirator	Standard/1.80 kg	TPK220
	Strong/2,70 kg	TPK330

Pack of 4 plastic Kirator retainers caps - Divergent.



NOT autoclavable – use cold steriliser. Maximum divergence of 44° between implants.

Example sequence









Ziacom®



ZM-Equator



ZM-Equator abutment

with applicator

Platf.	Height (H)	Reference
	1.00	ZMEX2001
	2.00	ZMEX2002
	3.00	ZMEX2003
	4.00	ZMEX2004
	5.00	ZMEX2005
	6.00	ZMEX2006
	1.00	ZMEX3401
	2.00	ZMEX3402
	3.00	ZMEX3403
	4.00	ZMEX3404
	5.00	ZMEX3405
	6.00	ZMEX3406
	1.00	ZMEX5001
	2.00	ZMEX5002
	3.00	ZMEX5003
	4.00	ZMEX5004

Golden surface treatment.



Includes ZM-Equator abutment with sterilisable polyoxymethylene plastic applicator (Tecaform AH-POM-C).

Related abutments

ZM-Equator impression transfer









System	Height (H)	Reference
ZM-Equator	6.50	TCRK3410



Pack of 4 units. DO NOT sterilise in an autoclave. Drillable.

ZM-Equator analogue





ZM-Equator

ZM-Equator processing kit



purple, transparent and pink.



System	Reference

ZM-Equator processing kit ZM8520 ZM-Equator processing kit consisting of: Titanium housing with black relined cap, spacer and three plastic caps in

Sterilise the metal coping in an autoclave. Plastic caps and the disc should be cold sterilised. See Instructions for Cleaning and Sterilising on the Ziacom® website.

System	Retention (kg)	Reference
	Light/1.20 kg	TZM100
ZM-Equator	Standard/1.80 kg	TZM200
	Strong/2.70 kg	TZM300

Pack of 4 plastic ZM-Equator retainer caps.



NOT autoclavable - use cold steriliser. Maximum divergence of 22° between implants.

ZM-Equator divergence processing kit





System	Reference
ZM-Equator processing kit	ZM8520D

ZM-Equator divergence processing kit comprising: Titanium housing with black relined cap, spacer and three plastic caps in purple, transparent and pink.

Sterilise the metal coping in an autoclave. Plastic caps and the disc should be cold sterilised. See Instructions for Cleaning and Sterilising on the Ziacom® website.

System	Retention (kg)	Reference
	Light/1.20 kg	TZM100
ZM-Equator	Standard/1.80 kg	TZM200
	Strong/2.70 kg	TZM300

Pack of 4 plastic ZM-Equator retainer caps - Divergent.



NOT autoclavable – use cold steriliser. Maximum divergence of 44° between implants.

Example sequence











ZM1

DIGITAL CAD-CAM

ZiaCam scanbody to implant







Platf.	Length (L)	Reference
	10.00	FNSYEX208T
	10.00	FNSYEX348T
	10.00	FNSYEX508T

Indicated for clinical and laboratory use.

All ZiaCam scanbodies to implant abutments include a screw Ref. LB102000 (NP)/LB103400 (RP/WP).

= = M1,60 M1,80

See the literature available at www.ziacom.com/biblioteca for more information on the use of zirconium restoration interfaces or the use of abutments in the "Prosthetic procedure" manual.



ZiaCam Ti-Base





Tx30 ZiaCam Ti-Base





Kirator. Toolbar abutment





Platf.	Height (H)	Reference
Universal	1.80	LOTB100

Golden surface treatment.





Rotatory

Non-rotatory


Platf.	Height (Hg/Ht)	Reference
	0.50/5.00	FRUEX201
	1.50/6.00	FRUEX202
	0.50/5.00	FRUEX341
	1.50/6.00	FRUEX342
	0.50/5.00	FRUEX501
	1.50/6.00	FRUEX502




F	Platf.	Height (Hg/Ht)	Reference
		0.50/5.00	FNUEX201
		1.50/6.00	FNUEX202
		0.50/5.00	FNUEX341
		1.50/6.00	FNUEX342
		0.50/5.00	FNUEX501
		1.50/6.00	FNUEX502

All ZiaCam Ti-Base abutments come with a special Kiran screw with surface treatment Ref. DSEI2010 (NP)/DSEI3410 (RP/WP).













Non-rotatory

Platf.	Height (Hg/Ht)	Reference
	0.50/6.00	FNUEX20TX1
	1.50/7.00	FNUEX20TX2
	0.50/6.00	FNUEX34TX1
	1.50/7.00	FNUEX34TX2
	0.50/6.00	FNUEX50TX1
	1.50/7.00	FNUEX50TX2
NO ROT	M1,80 M2,00	Grade 5 ELI Titanium















All Tx30 ZiaCam Ti-Base abutments come with a special Kiran Tx30 screw with surface treatment Ref. DSEI2010TX (NP)/DSEI3410TX (RP/WP).

Abutments Restorations using transepithelials




Restorations using transepithelials

■ Transepithelial abutments

- Allows the peri-implant tissue to form from the initial 8 weeks.
- · One abutment-one time, allows gingival adhesion to its surface as repeated detachments are not necessary.
- Avoids bone and soft tissue loss as there is no mechanical rupture of the peri-implant interface.
- The prosthetic working area is above the gingival level, making the soft tissue adhesive behaviour more predictable, maintaining a good seal.
- Less formation of micro-gaps at the implant-prosthesis junction.
- Increased crestal bone preservation.
- Prosthetic try-ins and definitive placement without anaesthesia.
- If the recommended torques are exceeded, the screw suffers the fracture at transepithelial level and not inside the implant.

Abutment heights


- · Greater abutment height means more marginal bone is preserved in cement-retained prostheses.
- Higher abutments (≥2mm) provide better soft tissue adaptation.
- Short abutments (< 2 mm) can compromise the soft tissues, resulting in more crestal bone loss.
- Marginal bone loss will differ depending on the clinical decision on the abutment height. Generally, prosthetic abutments ≥2mm will lead to better preservation of crestal bone.



7 30 Ziacom[®]



■ Basic | Example of usage sequence



ZM1 31 **Z**

Basic abutment





Platf.	Height (H)	Reference
	2.00	BASIC2002
	3.00	BASIC2003
	4.00	BASIC2004
	5.00	BASIC2005
	1.50	BASIC3401
	2.00	BASIC3402
	3.00	BASIC3403
	4.00	BASIC3404
	5.00	BASIC3405
	2.00	BASIC5002
	3.00	BASIC5003
	4.00	BASIC5004
	5.00	BASIC5005

Insertion key Ref. MABA100/MABA110.













Includes the Basic abutment with sterilisable polyoxymethylene applicator (Tecaform AH-POM-C). 18° cone angle. 36° angle between abutments.

Basic abutment with applicator

Basic healing abutment





System	Height (H)	Reference
Basic	5.00	BAHAEX34

Anodised —



















Basic impression abutment



System	Height (H)	Reference
Basic	8.00	BATC134


Anodised -



All Basic impression abutments come with a screw.

Basic analogue





System	Length (L)	Reference
Basic	13.00	BAIAEX34



Basic 3D analogue - Individual





System	Length (L)	Reference
Basic	13.00	BAIA348D



Basic 3D analogue - Pack





*Includes base screw Ref. DSIADI and lateral screw Ref. DSIADT for analogue connection.

Screws - 3D analogue







Туре	Length (L)	Reference	
Base screw (1)*	-	DSIADI	
Lateral screw (2)*	15.00	DSIADT	




*Pack of 4 units.




















ROT





ZM1 33 **Z**

DIGITAL CAD-CAM


ZiaCam scanbody to Basic abutment



Indicated for clinical and laboratory use.

All ZiaCam scanbodies to Basic abutments include a screw Ref. BDSEI3401.

ZiaCam Ti-Base to Basic



All Ti-Base ZiaCam to Basic abutments come with a special Kiran screw with surface treatment Ref. BDSEI3410.

ZiaCam Ti-Base Tx30 to Basic




All ZiaCam Ti-Base Tx30 to Basic abutments come with a special Kiran Tx30 screw with surface treatment Ref. BDSEI34TX.

Z 34 Ziacom®



Restorations using transepithelials

■ XDrive | Example of usage sequence



ZM1 35 **Z**

Abutments

XDrive straight abutment





Platf.	Height (H)	Reference
	1.00	XST103410
	2.00	XST103420
	3.00	XST103430
	4.00	XST103440
	5.00	XST103450

Insertion key Ref. MABA200/MABA210













Includes XDrive abutment with sterilisable polyoxymethylene applicator (Tecaform AH-POM-C).

21° cone angle. 42° angle between abutments.



XDrive abutment with applicator

XDrive 17° angled abutment







XDrive 30° angled abutment








latf.	Height (H)	Reference	Platf.	Height (H)	Reference
	2.00	XA2103417		3.00	XA3103430
	3.00	XA3103417		4.00	XA4103430
	4.00	XA4103417		5.00	XA5103430
	5.00	XA5103417	NO		Grade 5 ELI
O OT	(1,25mm) (1,80) (M2,00)	Grade 5 ELI Titanium	ROT	30 Ncm 1,25mm M1,80 M2,00	45° Titanium

All XDrive angled abutments come with a titanium positioner and screw.

XDrive healing abutment







XDrive impression abutment



Anodised




Includes screw.


XDrive analogue






XDrive 3D analogue - Pack








Screws - 3D analogue



*Pack of 4 units.











Kiran XDrive clinical screw



Special Kiran screw with surface treatment.

XDrive laboratory screw



NOT suitable for use as the final clinical screw.

Kiran Tx30 XDrive clinical screw





Reference

XST3410

For ZiaCam Ti-Base or metal structures




Kiran Tx30 special screw with surface treatment.

Length (L)

9.50

XDrive provisional abutment

XDrive provisional abutment










System




XDrive UCLA abutment





Machined base abutment XDrive

+ Castable abutment





Kirator XDrive abutment







37 🗷 ZM1

Abutments

DIGITAL CAD-CAM


ZiaCam scanbody to XDrive abutment



Indicated for clinical and laboratory use.


All ZiaCam scanbodies to XDrive abutments include a screw Ref. XLB103410.

ZiaCam XDrive Ti-Base



Includes special Kiran screw with surface treatment Ref.

ZiaCam Ti-Base Tx30 XDrive



Includes special Kiran Tx30 screw with surface treatment

■ Table of abutment torques

Element/Abutment	Instrument/Tool	Torque
Cover screws/Healing abutments	Hex screwdriver 1.25 mm	Manual
Impression abutment screws	Hex screwdriver 1.25 mm	Manual
Laboratory screws	Hex screwdriver 1.25 mm	Manual
Direct-to-implant clinical screws	Hex screwdriver 1.25 mm	30 Ncm
Kiran direct-to-implant clinical screws	Hex screwdriver 1.25 mm	30 Ncm
Nature abutments	Insertion keys: MANA100/MANA110/MANA120	30 Ncm
Clinical screws on Nature	Hex screwdriver 1.25 mm	30 Ncm
Kiran clinical screws on Nature	Hex screwdriver 1.25 mm	30 Ncm
Basic abutments	Insertion keys: MABA100/MABA110/MABA120	30 Ncm
XDrive abutments	Insertion keys: MABA200/MABA210/MABA220	30 Ncm
Clinical screws on Basic	Hex screwdriver 1.25 mm	25 Ncm
Kiran clinical screws on Basic	Hex screwdriver 1.25 mm	25 Ncm
Clinical screws on XDrive	Hex screwdriver 1.25 mm	20 Ncm
Kiran clinical screws on XDrive	Hex screwdriver 1.25 mm	20 Ncm
ZiaCam scanbody + screw	Hex screwdriver 1.25 mm	Manual
Kirator abutments	Insertion keys: LOSD01/LOSD02	30 Ncm
Tx30 abutment/screw (variable rotation)	Torx. screwdriver Tx30	30 Ncm

Exceeding the recommended tightening torque for screws and abutments compromises the prosthetic restoration and could damage the implant structure.

For immediate loading: DO NOT tighten manually, attach with the final torque. When using a screwdriver or adaptor for a contra-angle handpiece (CA), do not exceed a maximum speed of 25 rpm.

Surgical instruments



Surgical instruments

Universal Box surgical box



■ Universal Box contents available

5	
Empty BOX4	150U
Complete BOX45	50UC



Material: Radel.

Ensure boxes do not touch the walls of the autoclave to avoid damage.






■Surgical k	it contents	BOX450UC
		X
REF	Description	8
SID001M	Lance drill. Ø2.00 mm. Millimeter.	•
OSPD20M	Pilot Drill Ø1.60/2.00 mm. Millimeter.	•
OTD201C	Pilot Drill P1. Millimeter.	•
OTD301C	Pilot Drill P2. Millimeter.	•
OTD401C	Pilot Drill P3. Millimeter.	•
OTD203C	Final surgical drill. F1	•
OTD303C	Final surgical drill. F2	•
OTD403C	Final surgical drill. F3	•
OTD503C	Final surgical drill. F4	•
OTD603C	Final surgical drill. F5	•
CLD34	Crestal surgical drill. Ø 4.10 mm.	•
CLD50	Crestal surgical drill. Ø5.10 mm.	•
ZMPD160	Calibrated drill stop. 1. H6 mm.	•
ZMPD170	Calibrated drill stop. 1. H7 mm.	•
ZMPD185	Calibrated drill stop. 1. H8.5 mm.	•
ZMPD110	Calibrated drill stop. 1. H10 mm.	•
ZMPD115	Calibrated drill stop. 1. H11.5 mm.	•
ZMPD113	Calibrated drill stop. 1. H13 mm.	•
ZMPD114	Calibrated drill stop. 1. H14.5 mm.	•
ZMPD260	Calibrated drill stop. 2. H6 mm.	•
ZMPD270	Calibrated drill stop. 2. H7 mm.	•
ZMPD285	Calibrated drill stop. 2. H8.5 mm.	•
ZMPD210	Calibrated drill stop. 2. H10 mm.	•
ZMPD215	Calibrated drill stop. 2. H11.5 mm.	•
ZMPD213	Calibrated drill stop. 2. H13 mm.	•
ZMPD214	Calibrated drill stop. 2. H14.5 mm.	•
МТАРЗЗМС	Surgical tap. Ø3.30 mm. Millimeter.	•
MTAP36MC	Surgical tap. Ø3.60 mm. Millimeter.	•
MTAP40MC	Surgical tap. Ø4.00 mm. Millimeter.	•
MTAP44MC	Surgical tap. Ø4.40 mm. Millimeter.	•
MTAP48MC	Surgical tap. Ø4.80 mm Millimeter.	•
MUR101MT	Probe/Paralleling pin Initial. Millimeter.	•
MUR201MT	Probe/Paralleling pin P1. Millimeter.	•
MUR301MT	Probe/Paralleling pin P2. Millimeter.	•
MUR401MT	Probe/Paralleling pin P3. Millimeter.	•
DEXT10	Drill extender	•
MESD	Screwdriver tip, Ø 1.25 mm Long.	•
LMSD	1.25 mm surgical screwdriver. Long.	•
SMSD	1.25 mm surgical screwdriver. Short.	•
TORK50	Regulable torque wrench	•

The surgical kit box does not include the insertion keys or ZPlus block key. These must be ordered separately, according to the implant system used, Zinic® MTX / ZM10 / ZM1.

41 2 ZM1

Surgical instruments



2 42



STOPS

Calibrated drill stop





MPD160 MPD170 MPD185 MPD110
MPD185 MPD110
MPD110
MPD115
MPD113
MPD114
1PD260
1PD270
1PD285
1PD210
1PD215
1PD213
1PD214
MPD100
,

* Complete pack of 14 calibrated stops.



TAPS

Surgical tap. CA/Manual





Platf.	Diameter (Ø)	Reference
	3.30	MTAP33MC
	3.60	MTAP36MC
	4.00 *	MTAP40MC
	4.40 *	MTAP44MC
	4.80 *	MTAP48MC

Millimeter: 8.5/10/11.5/13/14.5

*Millimeter: 6/7/8.5/10/11.5/13/14.5




See surgical drilling protocol for more information on using tap.

ZM1

PROBES

Probe/Paralleling pin





	w			1.102
Platf.	Type	(W FWZ)	_	Reference
	Initial	1.60/2.00	27.00	MUR101MT
	Pilot 1	2.20/3.00	27.00	MUR201MT
	Pilot 2	2.60/3.50	27.00	MUR301MT
	Pilot 3	2.90/4.15	27.00	MUR401MT

Millimeter: 6/7/8.5/10/11.5/13/14.5



43 🗷

Surgical instruments

KEYS

ZPlus insertion key. Ratchet



Plati.	Length (L)	Reference
	3.10/Mini	XSMIN *
ZPlus	5.60/Short	TSMIN*
	10.60/Long	TLMIN*


● Hexagonal 2.4 mm / ■ Square 4x4 mm



*Ref. XSMIN/TSMIN/TLMIN are NOT included in the surgical

ZPlus insertion key. CA





Hexagonal 2.4 mm

ZPlus



*Ref. 01MMIN/02MMIN are NOT included in the surgical box.

Length (L) 15.90

23.90

ZM4 insertion key. Ratchet





Platf.	Length (L)	Reference
	15.00	SMEX20*
	15.00	SMEX34*
	15.00	SMEX50*

- Hexagonal NP 2.30 mm
- Hexagonal RP 2.70 mm
- Hexagonal WP 3.00 mm
- Square 4x4 mm



*Ref. SMEX20/SMEX34/SMEX50 are NOT included in the surgical box.

ZM4 insertion key. CA





Platf.	Length (L)	Reference
	7.50	MMEX20*
	7.50	MMEX34*
	7.50	MMEX50*

- Hexagonal NP 2.30 mm
- Hexagonal RP 2.70 mm
- Hexagonal WP 3.00 mm



*Ref. MMEX20/MMEX34/MMEX50 are NOT included in the surgical box.

Drill extender





Platf.	Length (L)	Reference
Universal	12.00	DEXT10





*Ref. 01MOHW are NOT included in the surgical box.

SCREWDRIVERS ZPlus block key Screwdriver tip. CA Surgical screwdriver. Manual Length (L) Platf. Reference Platf. Length (L) Reference Platf. Length (L) Reference 20.00/Short MESD01 2.80/Mini XSMSD Universal 90.00 01M0HW* **ZPlus** 25.00/Long MESD 9.50/Short SMSD Universal LMSD Hexagonal 1.25 mm 14.50/Long Hexagonal 2.4 mm 27.00/Extralong XLMSD

ZM1

Hexagonal 1.25 mm

RATCHET

Regulable torque wrench




Square 4x4 mm




45 🗷

Surgical instruments

Complementary instruments





Platf. Model Reference

ZM1 PRADIO60

Scales 1:1 and 1:1.25

Material: transparent acetate. Non-sterilisable material




2 46 Ziacom®

Prosthetic instruments



Prosthetic instruments

Prosthetic box



■ Contents of prosthetic boxes available

Contents	Reference
Empty	BOXPN
Basic	BOXPSN
Complete	BOXPCN



Material: Radel.

Ensure boxes do not touch the walls of the autoclave to avoid damage.



BOXPSN REF Description LOSD01 Kirator insertion key. MABA100 Basic insertion key. Short. MABA200 XDrive insertion key. Short. MADW10 Screwdriver handle. 4x4. SMSD1 Screwdriver tip. 1.25 mm. Short. LMSD1 Screwdriver tip. 1.25 mm. Long. XLMSD1 Screwdriver tip. 1.25 mm. Extra long. Screwdriver tip. 1.25 mm. Long. MESD Screwdriver tip. 1.25 mm. Short. MESD01 MESDTX Tx30 screwdriver tip. Long. LMSD1TX Tx30 screwdriver tip. Long. EDSZ20* ZPlus extractor screw. NP ZPlus extractor screw. RP/WP. EDSZ34*

Abutment extractor screw. RP

Abutment extractor screw. WP

Regulable torque wrench

EDSG34*

EDSG50*

TORK50

■ Contents of prosthetic boxes

^{*} Product not included in the ZM1 system.



KEYS

Kirator insertion key





System	Length (L)	Reference
Kirator	13.60/Ratchet/Manual	LOSD01
		LOSD02*

◆ Square 2.11 mm / ■ Square 4x4 mm



*Ref. LOSD02 is NOT included in the prosthetic box.

Basic insertion key. Ratchet





System	Length (L)	Reference
Basic	5.00/Short	MABA100
	13.00/Long	MABA110 *

Basic / Square 4x4 mm



*Ref. MABA110 is NOT included in the prosthetic box.

XDrive insertion key Ratchet





System	Length (L)	Reference
XDrive	6.00/Short	MABA200
	13.00/Long	MABA210*

OXDrive / Square 4x4 mm




*Ref. MABA210 is NOT included in the prosthetic box.

Nature insertion key. Ratchet







Nature / Square 4x4 mm



*Ref. MANA100/MANA110 are NOT included in the prosthetic

Nature insertion key. CA





*Ref. MABA210 is NOT included in the prosthetic box.








Basic insertion key. CA







Basic





*Ref. MABA210 is NOT included in the prosthetic box.

XDrive insertion key CA





MABA220*

XDrive **O** XDrive



*Ref. MABA220 is NOT included in the prosthetics box.

ZM1

7.00

49 🗷

Prosthetic instruments

SCREWDRIVERS


Screwdriver adapter handle



Square 4x4 mm



Screwdriver tip. Ratchet




Square 4x4 mm



Screwdriver tip. CA Length (L) Platf. Reference 20.00/Short MESD01 Universal **MESD** 25.00/Long



Tx30 screwdriver tip. CA



Do not exceed 30 Ncm as it could cause severe damage to the screwdriver and screw

* Ref. MESD01TX is NOT included in the prosthetics box.

Tx30 screwdriver tip. Ratchet





Do not exceed 30 Ncm as it could cause severe damage to

*Ref. SMSD1TX is NOT included in the prosthetic box.

Tx30 prosthetic screwdriver. Manual





Do not exceed 30 Ncm as it could cause severe damage to the screwdriver and screw


*Ref. SMSDTX/LMSDTX/XLMSDTX are NOT included in the prosthetics box.

EXTRACTOR SCREW

ZPlus extractor screw







Anodised NP RP/WP



* Product not included in the ZM1 system.

Galaxy/ZV2 abutment extractor screw



RATCHET

Regulable torque wrench





7 50 Ziacom®



Complementary instruments

CA to ratchet adaptor




Platf.	Length (L)	Reference
Universal	12.00	MC10Z

Square 4x4 mm



NOT included in the prosthetic box.

Extractor + Retainer inserter handle



NOT included in the prosthetic box.

Retention inserter



Platf.	Length (L)	Reference
Kirator	32.00	MBEI3602
ZM-Equator	32.00	MBEI3603



Kirator / ZM-Equator plastic cap insertion tool. NOT included in the prosthetic box.


Retentive joints instruments



Platf.	Dimensions	Reference
Universal	2x1	RREI0030

Pack of 10 units.

ZM1 51 **Z**



ZM1

Surgical protocols




Surgical protocol

Features of the ZM1 drilling system

■ Ziacom® drill system

Ziacom® implant system drills are made from stainless steel. A laser marking on the bur's shank identifies its inner and outer diameters and its length, while the horizontal laser marked bands on the active section corresponds to the different lengths of the implants (millimeter drills). The drill tip is 0.5 mm long and this is not included in the different laser-marked lengths.



■ Ziacom® final drills


Its use is essential and mandatory in order to achieve an ideal finish of the prepared implant bed for smooth, safe and precision insertion. In this way, overtorquing of the implant can be avoided while it is placed into its final position.

■ FINAL DRILL STOP

A stop, consisting of three blades (see red areas marked on image) has been incorporated into the design of the final drills, between the active area and the shank, to limit the penetration of the drill.

IMPORTANT

Take care not to drill beyond the stop, as this modifies the coronal anatomy of the surgical site.



■ ZIACOM® DRILLS EFFICIENCY GUARANTEE

Surgical drills for ZM1 implants from Ziacom® (cortical drills, lance drill, initial drill, pilot drills and final drills), have a lifetime of up to 60 uses. It is advisable to monitor the cutting status at all times, especially when reaching around 41 to 50 uses, since after 50 uses it is necessary to consider changing the drills before reaching 60 uses.

Bear in mind that, depending on the size of the implant, bone density and your surgical protocol, not all of the various drills will be used equally - it is recommended that you monitor the number of uses for each instrument.

40 **EFFICIENT USE** USE USES **50** 41 MONITORING CUTTING USES USES 51 60 CONSIDER CHANGING USES 60 LIFETIME - MAXIMUM USES

7 | 54




■ Ziacom® taps

Thread taps are available for contra-angle handpieces. The laser marking on the tap's shank identifies its diameter, while the horizontal laser marked bands on the active section corresponds to the different lengths.



■ Probe

Check the depth of the surgical site, especially when not using drill stops. To check the surgical bed axis, the paralleling pins are available in different diameters according to the drilling sequence.



■ Short and long insertion keys for ratchets and contra-angle handpieces

The insertion keys for contra-angle handpieces or ratchets have been designed for transporting implants from their No-Mount vial to the surgical site ready for insertion.





■ Drill stops

These are a surgical accessory that attach to drills and facilitate the work as they determine the depth of the osteotomy, providing extra assurance when preparing the surgical site.




ZM1 55 **Z**

Surgical protocol

Features of the ZM1 drilling system

■ Internal view of the Universal Box surgical box



Recommendations on the maximum implant insertion torque



The recommended insertion torque ranges between 35 and 50 Ncm on a case-by-case basis.

To avoid deforming the driver and/or implant connection, insertions performed with a contra-angle handpiece (CA) must respect the recommended maximum rpm (25 rpm) and maximum torque (50 Ncm).

If the implant cannot be fully inserted using the recommended maximum torque, withdraw the implant, repeat the drilling and then re-insert it.

Check the final insertion torque with the adjustable dynamometric ratchet Ref. TORK50 or a contra-angle handpiece.

Exceeding the maximum torque (50 Ncm) when inserting the implant may result in:


- · Irreversible deformations in the implant's internal connection.
- · Irreversible deformations in the implant insertion instruments.
- Difficulty or impossibility in dismounting the instrument/implant assembly.

Z 56 Ziacom[®]



■ ZM1 implant

It is important to note that the drilling protocol for ZM1 implants using stepped drill bits varies significantly based on the implant diameter and the type of bone at the surgical site, and therefore it is important to pay special attention to these two aspects.



High-density drilling protocol steps (D1 - D2*)

PRELIMINARY STEP | Opening the gum

Make an incision and raise the flap.



STEP 1 | Lance drill



Start the implant site drilling sequence using the Lance Drill Ref. SID001M. Be aware of the laser marking on the drill to indicate the length, or use drill stop Ref. ZMPD115. Control the direction and angle of drilling by applying intermittent pressure vertically, taking care not to exert too much pressure on the bone. If necessary, use drill extender Ref. DEXT10.





STEP 2 | Initial drill



Continue the drilling sequence using Initial Drill Ref. OSPD20M until the total length of the chosen implant is reached. Be aware of the laser marking on the drill that indicates the length, or use the drill stop Ref. ZMPD115. Monitor the direction and inclination of the drilling, exerting pressure intermittently, always in a vertical direction, taking care not to generate excessive pressure on the bone. If necessary, use drill extender Ref. DEXT10.





STEP 3 | Initial Drill Probe/Paralleling pin



Check the depth of the surgical site and the insertion axis by inserting the Initial Drill Probe/Paralleling pin Ref. MUR101MT. Repeat this step as many times as necessary during the surgery.

ZM1 57

Surgical protocol

STEP 4 | Pilot drill 1



Continue the drilling sequence using Pilot Drill 1 Ref. OTD201C, until the full length of the chosen implant is reached. Be aware of the laser marking on the drill to indicate the length, or use drill stop Ref. ZMPD115. Control the direction and angle of drilling by applying intermittent pressure vertically, taking care not to exert too much pressure on the bone. If necessary, use drill extender Ref. DEXT10.





NOTE

Once this step has been completed, to fit an implant with diameter:

- Ø3.30 mm > Final drill 1 Ref. OTD203C + Tap MTAP33MC
- Ø3.60 mm > Final drill 2 Ref. OTD303C + Tap MTAP36MC

STEP 6 | Pilot drill 2



Continue the drilling sequence using Pilot Drill 2 Ref. OTD301C, until the full length of the chosen implant is reached. Be aware of the laser marking on the drill to indicate the length, or use drill stop Ref. ZMPD115. Control the direction and angle of drilling by applying intermittent pressure vertically, taking care not to exert too much pressure on the bone. If necessary, use drill extender Ref. DEXT10.





NOTE

Once this step has been completed, to fit an implant with diameter:

- Ø4.00 mm > Final Drill 3 Ref. OTD403C) + Tap MTAP40MC
- Ø4.40 mm > Final Drill 4 Ref. OTD503C) + Tap MTAP44MC

STEP 8 | Final Drill 3



Continue the drilling sequence using Pilot Drill 3 Ref. OTD403C up to the length corresponding to the cortical bone thickness, according to individual clinical case. Control the direction and angle of drilling by applying intermittent pressure vertically, taking care not to exert too much pressure on the bone. If necessary, use drill extender Ref. DEXT10.





STEP 5 | Pilot Probe/Paralleling pin 1



Check the depth of the surgical site and the insertion axis by inserting the Pilot Probe/Paralleling pin 1 Ref. MUR201MT. Repeat this step as many times as necessary during the surgery.

STEP 7 | Pilot Probe/Paralleling pin 2



Check the depth of the surgical site and the insertion axis by inserting the Pilot Probe/Paralleling pin 2 Ref. MUR301MT. Repeat this step as many times as necessary during the surgery.

STEP 9 | Surgical tap Ø4.00



Place the Ø4.00 mm surgical tap, Ref. MTAP40MC in the surgical site. Apply firm pressure and start to turn slowly. Once threads engage, continue to screw the tap in without pressure to the planned depth. If excessive resistance is met, turn 90° anti-clockwise after each complete turn. To remove the tap, turn it anti-clockwise. While using the tap, it is recommended that you pass it along the entire length of the implant.





■ Important notes: Type D2* Bone Density

In the case of type D2 bone density, the surgical drilling protocol indicated for type D1 bone density should be followed, leaving out the use of the Surgical Tap on any of the implant diameters. Nevertheless, it is up to the discretion of the professional to decide on full or partial use the Surgical Tap, based on their clinical experience and the identification of the density of the existing bone at the site. This is particularly relevant in cases where the bone density varies significantly along the length of the osteotomy for the implant.


ZM1 59 **Z**

Surgical protocol

Low-density drilling protocol steps (D3 - D4**)

PRELIMINARY STEP | Opening the gum

Make an incision and raise the flap



STEP 1 | Lance Drill



Start the implant site drilling sequence using the Lance Drill Ref. SID001M. Be aware of the laser marking on the drill to indicate the length, or use drill stop Ref. ZMPD115. Control the direction and angle of drilling by applying intermittent pressure vertically, taking care not to exert too much pressure on the bone. If necessary, use drill extender Ref. DEXT10.





STEP 2 | Initial drill



Continue the drilling sequence using Initial Drill Ref. OSPD20M until the total length of the chosen implant is reached. Be aware of the laser marking on the drill that indicates the length, or use the drill stop Ref. ZMPD115. Monitor the direction and inclination of the drilling, exerting pressure intermittently, always in a vertical direction, taking care not to generate excessive pressure on the bone. If necessary, use drill extender Ref. DEXT10.





Once this step has been completed, to fit an implant with diameter:

- Ø3.30 mm > Final Drill 1 Ref. OTD203C)
 Ø3.30 mm > Final Drill 2 Ref. OTD303C)

STEP 3 | Initial Drill Probe/Paralleling pin



Check the depth of the surgical site and the insertion axis by inserting the Initial Drill Probe/Paralleling pin Ref. MUR101MT. Repeat this step as many times as necessary during the surgery.

STEP 4 | Pilot drill 1



Continue the drilling sequence using Pilot Drill 1Ref. OTD201C, until the full length of the chosen implant is reached. Be aware of the laser marking on the drill to indicate the length, or use drill stop Ref. ZMPD115. Control the direction and angle of drilling by applying intermittent pressure vertically, taking care not to exert too much pressure on the bone. If necessary, use drill extender Ref. DEXT10.





Once this step has been completed, to fit an implant with diameter:

- Ø3.30 mm > Final Drill 3 Ref. OTD403C)
- Ø4.40 mm > Final Drill 4 Ref. OTD503C)

STEP 5 | Pilot Probe/Paralleling pin 1



Check the depth of the surgical site and the insertion axis by inserting the Pilot Probe/Paralleling pin 1 Ref. MUR201MT. Repeat this step as many times as necessary during the surgery.



STEP 7 | Final Drill 3



Continue the drilling sequence using Pilot Drill 3 Ref. OTD403C up to the length corresponding to the cortical bone thickness, according to individual clinical case. Control the direction and angle of drilling by applying intermittent pressure vertically, taking care not to exert too much pressure on the bone. If necessary, use drill extender Ref. DEXT10.





■ Important notes: Type D4** Bone Density

In the case of type D4 bone density, the surgical drilling protocol indicated for type D3 bone density should be followed, leaving out the use of the last Final Drill for each of the implant diameters. Nevertheless, it is up to the discretion of the professional to decide to use the last Final Drill fully or partially, based on their clinical experience and the identification of the density of the existing bone at the site. This is particularly relevant in cases where the bone density varies significantly along the length of the osteotomy for the implant.


Surgical protocol

Implant placement with ZPlus Mount | Titansure

ZPlus Mount

Surface treatment

Titansure



STEP 1 | Unpacking the implant

- 11 Press the word "PRESS" and tear open the hox
- Remove the top of the carton and take out the blister pack.
- (13) Carefully remove the seal from the blister pack.
- 14 Turn the vial containing the implant out onto a sterile cloth in the operating area.
- 15 Remember to remove the label from the implant and to stick it onto the patient's implant card and clinical records to ensure that the product is traceable.




STEP 2 | Choosing the right insertion instrument

Based on the specific clinical situation and access to the surgical site, one of three different instruments can be selected to insert the implant:

- (A) Contra-angle: use the ZPlus insertion key. CA driver of the desired length Ref. 01MMIN / 02MMIN and insert it into the contra-angle.
- (3) Torque wrench Ref. TORK50: use the ZPlus insertion key. Ratchet/Manual of the desired length Ref. XSMIN / TSMIN / TLMIN and insert it into the ratchet set to function "IN", which is identified with an arrow.
- © Screwdriver handle 4x4 Ref. MADW10: use the ZPlus insertion key. Ratchet/ Manual of the desired length Ref. XSMIN / TSMIN / TLMIN and insert it into the screwdriver handle.

STEP 3 | Removing the implant from its vial


Hold the vial containing the implant in one hand and insert the selected ZPlus insertion key with the other hand. Remove the implant-mount assembly by lifting it vertically out of the vial.





ZM1 implant insertion with ZPlus Mount

STEP 4 | Inserting the implant




Insert the implant into the surgical site, controlling both the direction and angle of the implant. When inserting the implant with a contra-angle, use a maximum speed of 25 rpm. The recommended insertion torque ranges from 35 to 50 Ncm, according to each case, and is not limited to a single torque.

If resistance is met during insertion, turn the implant slightly anti-clockwise and then continue to insert after waiting a few seconds. Repeat this process as many times as necessary.

The Ziacom® surgical protocol establishes crestal positioning of the implant platform.

The ZPlus has 3 flat sides. After inserting the implant, make sure that one of these flat sides faces the vestibular direction.



Use direct-to-implant insertion keys, Ref. SMEX20/SMEX34/SMEX50 for Torque Wrench/Manual and MMEX20/MEX34/MMEX50 for CA, to adjust the final position of the implant.

STEP 5A | Extracting the ZPlus Mount

Lock the ZPlus mount using block key Ref. 01MOHW and remove the screw using manual surgical screwdriver Ref. SMSD / LMSD. After removing the screw, remove the ZPlus by hand.




REF. SMSD/LMSD with 01MOHW

IMPORTANT

If the ZPlus gets caught on the implant, use the extractor screw to remove it: with NP platform, use Ref. EDSZ20 and with RP/WP platform, use Ref. EDSZ34.

STEP 5B | Extracting the ZPlus Mount

In order to prevent the ZPlus mount from warping or cold welding with the implant, the point of insertion at which the mount should be extracted will depend on the type of bone.





IMPORTANT

The maximum insertion torque for the dental implants is 50 Ncm. Exceeding the maximum insertion torque for the implants may cause severe damage to the dental implant, its connection, the Mount and the clinical screw included. Check the specifications in the surgical protocol for removal of the Mount, according to the type of implant connection and the bone type.

STEP 6 | Crestal placement of the implant

The Ziacom® ZM1 implant platform should be placed at bone ridge level.



RECOMMENDED ridge position

ZM1 63 **Z**

Surgical protocol

■ Soft tissue conditioning

STEP 1 | Placing the cover screw



Remove the cover screw anti-clockwise using manual surgical screwdriver Ref. SMSD / LMSD. Move the cover screw towards the implant while taking care not to drop it and cause its accidental ingestion. Insert the screw into the implant until it locks, applying manual torque in a clockwise direction. Placement of the cover screw during the first surgical phase requires that, after the osseointegration period, the second surgical phase should be performed or the implant should be exposed to fit the chosen abutment.

Based on each individual case, you can choose not to place a cover screw but instead to directly attach a healing abutment.




STEP 2 | Closing the soft tissue

Close and suture the soft tissue, carefully lining up the flaps.



STEP 3 | Exposing and extracting the cover screw



STEP 4 | Placing the healing abutment



Insert the chosen healing abutment using manual surgical screwdriver Ref. SMSD or LMSD.

The choice of healing abutment will depend on each individual case. It should match the implant platform and also the height of the gingival tissue in order to prevent occlusion of the abutment. If the abutment is too tall, it may subject the implant to premature loading, compromising the osseointegration process.



2 64 Ziacom®



Bone types

Misch classification (1988)



TYPE **D1** BONE

- Dense cortical and dense trabecular
 bone
- > 1250 HU



TYPE **D2** BONE

- Porous cortical and dense trabecular bone.
- 850 1250 HU



TYPE **D3** BONE

- Porous cortical and fine trabecular bone.
- 350 850 HU



TYPE **D4** BONE

- Thin crestal cortical and fine trabecular bone.
- 150 350 HU


HU = Hounsfield Units

IMPORTANT

In order to simplify the surgical drilling protocols, we have created quick drilling guides, in which the criteria for bone types are amalgamated, with D1-D2 treated as "High-Density" bone, and D3-D4 bone types as "Low-Density" bone.

■ Handling of cover screw

Place the cover screw in the screwdriver. Move the cover screw towards the implant while taking care not to drop it and cause its accidental ingestion. Insert it into the implant applying manual torque in a clockwise direction.



Considerations for temporisation and immediate loading

Immediate temporisation and immediate loading are procedures that involve the placement of the prosthesis within 72 hours after implant surgery. The fundamental difference between these procedures is whether or not the prosthesis will have a functional load.

Adequate primary stability of the implant at the time of insertion is crucial to consider placing a provisional or immediately loaded prosthesis. This stability can be objectively measured by the insertion torque, which must be equal to or greater than 40-45 Ncm or by analysing the resonance frequency (ISQ value), which should be greater than or equal to 70.

■ IMMEDIATE TEMPORISATION

Immediate temporisation involves thorough monitoring of occlusion, both in central (closed) position, and during lateral or dynamic movements that occur during mastication. By freeing the provisional from any contact in these situations, the transfer of forces to the implant is prevented.

The main objectives of immediate temporisation are:

- Immediate closure of edentulous spaces in aesthetic areas.
- · Guided regeneration of the gingival emergence profile due to the presence of the provisional crown or bridge.

■ IMMEDIATE LOADING

The principle of immediate loading involves the controlled transfer of contact from the moment of placement of the restoration while the restoration is in occlusion; therefore we distinguish between:

- · Progressive immediate loading, using an acrylic provisional restoration as the initial restoration (released in dynamic occlusion).
- Definitive immediate loading, with rigid material and active occlusion from day one.

Both processes involve risks to the success of the osseointegration of the implant, so it is up to the practitioner, based on clinical experience and the case in question, whether or not to place an immediate provisional restoration and/or immediate loading.

ZM1 65 **Z**

Surgical protocol

Restorations using transepithelials

■ Transepithelial abutments

- Allows the peri-implant tissue to form from the initial 8 weeks.
- · One abutment-one time, allows gingival adhesion to its surface as repeated detachments are not necessary.
- Avoids bone and soft tissue loss as there is no mechanical rupture of the peri-implant interface.
- The prosthetic working area is above the gingival level, making the soft tissue adhesive behaviour more predictable, maintaining a good seal.
- · Less formation of micro-gaps at the implant-prosthesis junction.
- · Increased crestal bone preservation.
- Prosthetic try-ins and definitive placement without anaesthesia.
- If the recommended torques are exceeded, the screw suffers the fracture at transepithelial level and not inside the implant.

Abutment heights

- · Greater abutment height means more marginal bone is preserved in cement-retained prostheses.
- Higher abutments (≥2mm) provide better soft tissue adaptation.
- Short abutments (< 2 mm) can compromise the soft tissues, resulting in more crestal bone loss.
- Marginal bone loss will differ depending on the clinical decision on the abutment height. Generally, prosthetic abutments ≥2mm will lead to better preservation of crestal bone.




Simplified surgical protocol

These surgical guides have been designed with a simplified surgical protocol to perform simple and efficient drilling of the surgical site.

ZPlus - Drilling Protocol




Detailed speeds are recommended



ZM1 67 **Z**

Surgical protocol

 $ZM1 \not 03.60$ (Example of preparation of the implant bed with ZM1 implant Ø3.60x11.5)






2 68




 $ZM1 \not 04.00$ (Example of preparation of the implant bed with ZM1 implant $\emptyset 4.00$ x11.5)



ZM1 69 **Z**


Simplified surgical protocol

 Z_{M1} Ø4,40 (Example of preparation of the implant bed with ZM1 implant Ø4.40x11.5)





 Zm_1 Ø4,80 (Example of preparation of the implant bed with ZM1 implant Ø4.00x11.5)





Simplified surgical protocol

General recommendations

■ Points to consider during the procedure

- Surgical drills must be inserted into the contra-angle handpiece with the motor stopped, ensuring that they are seated and rotate properly before starting drilling. Treat drills with the utmost care; the slightest damage to the tips could compromise their effective operation.
- Damaged instruments must be disposed of according to local regulations.
- Implantologists should keep one of the identification labels supplied with the product in the patient's records so that the product can be traced correctly.
- Each instrument must only be used for the specific use recommended by the manufacturer.

Always consult the surgical and prosthetic protocols published in this catalogue, as well as the other documents available in the "Reference literature" section of our website **www.ziacom.com/biblioteca** which explain the procedures, protocols and instructions for use before using the Ziacom ZM1 system.



72 Ziacom[®]

Cleaning, disinfection and sterilisation



Cleaning, disinfection and sterilisation

The protocols described in this section must only be carried out by personnel qualified to clean, disinfect and sterilise the dental materials specified herein.

Cleaning and disinfection instructions

Applicable for surgical and prosthetic instruments and boxes.

Disassembly

- 1. Disassemble* the instruments that need to be cleaned and disinfected, such as manual ratchets, drills or drill stops.
- 2. Remove all the different components from the surgical or prosthetic kit box for correct cleaning.

Cleaning and disinfection

For disinfection of instruments and surgical kit boxes:

- 1. Submerge the instruments in a detergent/disinfectant solution** suitable for dental instruments to help eliminate any adhered biological residues. If an ultrasound bath is available***, confirm that the detergent/disinfectant solution is indicated for use with this type of equipment.
- 2. Manually remove any biological residues with a non-metallic brush and pH-neutral detergent.
- 3. Rinse with copious water.
- 4. When cleaning surgical and prosthetic kit boxes, always use a pH-neutral detergent and non-abrasive tools to avoid damaging the surface of the boxes.
- 5. Dry the materials with disposable, lint-free, cellulose cloths or compressed air.

For disinfection of plastic caps and the protective disk:

- 1. Submerge for 10 minutes in a neat benzalkonium chloride solution.
- 2. Rinse with distilled water.
- 3. Dry the caps and disk prior to use.

Inspection

- 1. Check that the instruments are perfectly clean; if not, repeat the cleaning and disinfection steps.
- 2. Discard any instruments with imperfections and replace them before the next surgery.
- 3. Check that the instruments and surgical and prosthetic kit boxes are perfectly dry before reassembling the parts and proceeding with sterilisation.
 - * See the assembly and disassembly manuals at www.ziacom.com/biblioteca
 - ** Follow the instructions from the disinfectant's manufacturer to determine the correct concentrations and times.
 - ** Follow the instructions from the ultrasound bath's manufacturer to determine the correct temperature, concentration and times.

Sterilisation instructions for steam autoclaves

Applicable to orthodontic implants, abutments, kit, surgical and prosthetic boxes, pins, fixing screws and mesh membranes.

- 1. Place the material in individual sterilisation pouches and seal the pouches. For joint sterilisation, place the instruments in their surgical kit box, place the box in a sterilisation pouch and seal the pouch.
- 2. Place the pouches to be sterilised in the autoclave.
- 3. Sterilise in a steam autoclave at 134°C/273°F (max. 137°C/276°F) for 4 min (minimum) at 2 atm. Dynamometric torque wrenches must be sterilised in 3 vacuum cycles at 132°C/270°F for at least ≥ 4 minutes and vacuum dried for at least 20 minutes.

For the United States only: The validated and recommended sterilisation cycle for the US must be performed in a steam autoclave at 132°C/270°F for at least 15 minutes with a drying time of at least 15–30 minutes.

IMPORTANT

Make sure the drying stage is allowed to run to completion, otherwise the products may be damp.

Check the sterilisation equipment if the materials or sterilisation pouches are damp at the end of the sterilisation cycle.

Perform the necessary maintenance actions on the autoclave according to the established periodicity and following the manufacturer's instructions.

Z



Storage of Ziacom® products

- Store the products in their original packaging in a clean, dry place until they are to be used.
- · After sterilisation, keep the products in the sealed sterilisation pouches in a clean, dry location.
- Never exceed the use by date indicated by the manufacturer of the sterilisation pouches.
- Always follow the instructions of the manufacturer of the sterilisation pouches.

General recommendations

- Never use damaged or dirty material; never reuse single-use products. The user is responsible for following the instructions described in this document correctly.
- · Pay attention to piercing or sharp elements. Gloves should be worn when cleaning the materials to avoid accidents during handling.
- Follow the safety instructions indicated by the manufacturer of the disinfectant.
- The product's sterility cannot be guaranteed if the sterilisation pouch is open, damaged or damp.
- Respect all stages of the sterilisation process. If the materials or sterilisation pouches contain traces of water or moisture, check the autoclave and repeat the sterilisation.
- Orthodontic abutments and implants are supplied UNSTERILISED and must always be sterilised before use.
- Instruments and surgical and prosthetic kit boxes are supplied UNSTERILISED and must always be sterilised before use and cleaned and disinfected after use.
- Sterilisation, cleaning and disinfection processes gradually deteriorate the instruments. Inspect the instruments thoroughly to detect any signs of deterioration.
- Avoid contact between products made from different materials (steel, titanium, etc.) during the cleaning, disinfection and sterilisation processes.
- Ziacom Medical SL recommends these instructions are implemented for the correct maintenance and safety of their products; accordingly, the company refuses any liability for any damage to the products that could arise if the user applies alternative cleaning, disinfection and sterilisation procedures.

See the latest version of the cleaning, disfection and sterilisation instructions at www.ziacom.com/biblioteca



ZM1 75 Z





See the updated general conditions of sale at www.ziacom.com.

Check the availability of each product in your country.

All rights reserved. No part of this document may be reproduced or stored in any medium or reproduction system, nor transmitted in any way or under any concept, electronically, mechanically, in photocopies, recording or any other mean not considered here without the permission of holder of the copyright, editing and printing. Ziacom® is a registered trademark of Ziacom Medical SL.

See the latest version of the catalogues at www.ziacom.com.

EN|ENGLISH





